


## Types of Auto-injectors

Adrenaline autoinjector	<b>EpiPen</b> 	<b>Emerade</b> 	<b>Jext</b> 
Website	<a href="#">EpiPen® &amp; EpiPen® Jr: Official UK Site</a>	<a href="#">Auto injector for anaphylaxis   Emerade (emerade-bausch.co.uk)</a>  Currently withdrawn from UK	<a href="#">Welcome to Jext</a>
Patient Information leaflet	<a href="#">pil.4289.pdf (medicines.org.uk)</a>	<a href="#">pil.5278.pdf (medicines.org.uk)</a>	<a href="#">JEXT (medicines.org.uk)</a>
Dosage (British National Formulary)	<a href="#">Adrenaline/epinephrine   Drugs   BNF   NICE</a>	<a href="#">Adrenaline/epinephrine   Drugs   BNF   NICE</a>	<a href="#">Adrenaline/epinephrine   Drugs   BNF   NICE</a>
Needle length	150mcg device: 13mm 300mcg device: 15mm	300mcg device: 25mm 500mcg device: 25mm	150mcg device: 13mm 300mcg device: 15mm
Demonstration video for adults	<a href="#">How To Use: Administering Your EpiPen®   EpiPen®</a>	<a href="#">How to use Emerade   Emerade (emerade-bausch.co.uk)</a>	<a href="#">Jext® Adults</a>
Demonstration video for children	<a href="#">How To Use: Administering Your EpiPen®   EpiPen®</a>	<a href="#">How to use Emerade   Emerade (emerade-bausch.co.uk)</a>	<a href="#">Jext® Kids</a>  <a href="#">Jext® Teens</a>
Expiry alert service	<a href="#">EpiPen® Expiry Date Alerts   EpiPen®</a>	<a href="#">Patient   Emerade (emerade-bausch.co.uk)</a>	<a href="#">Jext® Adults</a>  <a href="#">Jext® Teens</a>  <a href="#">Jext® Kids</a>
Trainer pens available	<a href="#">EpiPen® Training Pen   EpiPen®</a>	<a href="#">Patient Resources   Emerade (emerade-bausch.co.uk)</a>	<a href="#">Order Trainer Pen (jext.co.uk)</a>